

SECTION 3 (3.01-3.10)

Small Fragment



Screws

VO.03.1
3.5mm Cortical Screws

Drill Bit : 2.5mm



Length (mm)	Stainless Steel	Titanium
10	VO.03.1.SS01	VO.03.1.TI01
12	VO.03.1.SS02	VO.03.1.TI02
14	VO.03.1.SS03	VO.03.1.TI03
16	VO.03.1.SS04	VO.03.1.TI04
18	VO.03.1.SS05	VO.03.1.TI05
20	VO.03.1.SS06	VO.03.1.TI06
22	VO.03.1.SS07	VO.03.1.TI07
24	VO.03.1.SS08	VO.03.1.TI08
26	VO.03.1.SS09	VO.03.1.TI09
28	VO.03.1.SS10	VO.03.1.TI10
30	VO.03.1.SS11	VO.03.1.TI11
32	VO.03.1.SS12	VO.03.1.TI12
34	VO.03.1.SS13	VO.03.1.TI13
36	VO.03.1.SS14	VO.03.1.TI14
38	VO.03.1.SS15	VO.03.1.TI15
40	VO.03.1.SS16	VO.03.1.TI16
42	VO.03.1.SS17	VO.03.1.TI17
44	VO.03.1.SS18	VO.03.1.TI18
46	VO.03.1.SS19	VO.03.1.TI19
48	VO.03.1.SS20	VO.03.1.TI20
50	VO.03.1.SS21	VO.03.1.TI21

VO.03.3
3.5mm Cancellous Screws

Drill Bit : 2.5mm



Length (mm)	Stainless Steel	Titanium
10	VO.03.3.SS01	VO.03.3.TI01
12	VO.03.3.SS02	VO.03.3.TI02
14	VO.03.3.SS03	VO.03.3.TI03
16	VO.03.3.SS04	VO.03.3.TI04
18	VO.03.3.SS05	VO.03.3.TI05
20	VO.03.3.SS06	VO.03.3.TI06
22	VO.03.3.SS07	VO.03.3.TI07
24	VO.03.3.SS08	VO.03.3.TI08
26	VO.03.3.SS09	VO.03.3.TI09
28	VO.03.3.SS10	VO.03.3.TI10
30	VO.03.3.SS11	VO.03.3.TI11
35	VO.03.3.SS12	VO.03.3.TI12
40	VO.03.3.SS13	VO.03.3.TI13
45	VO.03.3.SS14	VO.03.3.TI14
50	VO.03.3.SS15	VO.03.3.TI15
55	VO.03.3.SS16	VO.03.3.TI16
60	VO.03.3.SS17	VO.03.3.TI17

VO.03.2
3.5mm Cortical Screws, Self Tapping

Drill Bit : 2.5mm



Length (mm)	Stainless Steel	Titanium
10	VO.03.2.SS01	VO.03.2.TI01
12	VO.03.2.SS02	VO.03.2.TI02
14	VO.03.2.SS03	VO.03.2.TI03
16	VO.03.2.SS04	VO.03.2.TI04
18	VO.03.2.SS05	VO.03.2.TI05
20	VO.03.2.SS06	VO.03.2.TI06
22	VO.03.2.SS07	VO.03.2.TI07
24	VO.03.2.SS08	VO.03.2.TI08
26	VO.03.2.SS09	VO.03.2.TI09
28	VO.03.2.SS10	VO.03.2.TI10
30	VO.03.2.SS11	VO.03.2.TI11
32	VO.03.2.SS12	VO.03.2.TI12
34	VO.03.2.SS13	VO.03.2.TI13
36	VO.03.2.SS14	VO.03.2.TI14
38	VO.03.2.SS15	VO.03.2.TI15
40	VO.03.2.SS16	VO.03.2.TI16
42	VO.03.2.SS17	VO.03.2.TI17
44	VO.03.2.SS18	VO.03.2.TI18
46	VO.03.2.SS19	VO.03.2.TI19
48	VO.03.2.SS20	VO.03.2.TI20
50	VO.03.2.SS21	VO.03.2.TI21

VO.03.4
4.0mm Cancellous Screw, Short Thread

Drill Bit : 2.5mm



Length (mm)	Stainless Steel	Titanium
10	VO.03.4.SS01	VO.03.4.TI01
12	VO.03.4.SS02	VO.03.4.TI02
14	VO.03.4.SS03	VO.03.4.TI03
16	VO.03.4.SS04	VO.03.4.TI04
18	VO.03.4.SS05	VO.03.4.TI05
20	VO.03.4.SS06	VO.03.4.TI06
22	VO.03.4.SS07	VO.03.4.TI07
24	VO.03.4.SS08	VO.03.4.TI08
26	VO.03.4.SS09	VO.03.4.TI09
28	VO.03.4.SS10	VO.03.4.TI10
30	VO.03.4.SS11	VO.03.4.TI11
35	VO.03.4.SS12	VO.03.4.TI12
40	VO.03.4.SS13	VO.03.4.TI13
45	VO.03.4.SS14	VO.03.4.TI14
50	VO.03.4.SS15	VO.03.4.TI15
55	VO.03.4.SS16	VO.03.4.TI16
60	VO.03.4.SS17	VO.03.4.TI17

Screws

VO.03.5
4.0mm Cancellous Screw, Full Thread

Drill Bit : 2.5mm



Length (mm)	Stainless Steel	Titanium
10	VO.03.5.SS01	VO.03.5.TI01
12	VO.03.5.SS02	VO.03.5.TI02
14	VO.03.5.SS03	VO.03.5.TI03
16	VO.03.5.SS04	VO.03.5.TI04
18	VO.03.5.SS05	VO.03.5.TI05
20	VO.03.5.SS06	VO.03.5.TI06
22	VO.03.5.SS07	VO.03.5.TI07
24	VO.03.5.SS08	VO.03.5.TI08
26	VO.03.5.SS09	VO.03.5.TI09
28	VO.03.5.SS10	VO.03.5.TI10
30	VO.03.5.SS11	VO.03.5.TI11
35	VO.03.5.SS12	VO.03.5.TI12
40	VO.03.5.SS13	VO.03.5.TI13
45	VO.03.5.SS14	VO.03.5.TI14
50	VO.03.5.SS15	VO.03.5.TI15
55	VO.03.5.SS16	VO.03.5.TI16
60	VO.03.5.SS17	VO.03.5.TI17

VO.03.6
3.5mm Malleolar Screw

Drill Bit : 2.5mm



Length (mm)	Stainless Steel	Titanium
25	VO.03.6.SS01	VO.03.6.TI01
30	VO.03.6.SS02	VO.03.6.TI02
35	VO.03.6.SS03	VO.03.6.TI03
40	VO.03.6.SS04	VO.03.6.TI04
45	VO.03.6.SS05	VO.03.6.TI05
50	VO.03.6.SS06	VO.03.6.TI06
55	VO.03.6.SS07	VO.03.6.TI07
60	VO.03.6.SS08	VO.03.6.TI08
65	VO.03.6.SS09	VO.03.6.TI09
70	VO.03.6.SS10	VO.03.6.TI10

VO.03.7
Washer for Small Screws



	Stainless Steel	Titanium
Code	VO.03.7.SS01	VO.03.7.TI01

SECTION 3

SMALL FRAGMENT

Plates

VO.03.8 3.5mm One-third Tubular Plate

Used as buttress plate for Tubular, Metatarsals and Metacarpals.
Fixation with 3.5 & 4.0mm Screws. Compression.



Holes	Stainless Steel	Titanium
2	VO.03.8.SS01	VO.03.8.TI01
3	VO.03.8.SS02	VO.03.8.TI02
4	VO.03.8.SS03	VO.03.8.TI03
5	VO.03.8.SS04	VO.03.8.TI04
6	VO.03.8.SS05	VO.03.8.TI05
7	VO.03.8.SS06	VO.03.8.TI06
8	VO.03.8.SS07	VO.03.8.TI07
9	VO.03.8.SS08	VO.03.8.TI08
10	VO.03.8.SS09	VO.03.8.TI09
11	VO.03.8.SS10	VO.03.8.TI10
12	VO.03.8.SS11	VO.03.8.TI11

VO.03.10 3.5mm One-third Tubular Plate with Collar

Used as buttress plate for Tubular, Metatarsals and Metacarpals.
Fixation with 3.5 & 4.0mm Screws. Compression.



Holes	Stainless Steel	Titanium
2	VO.03.10.SS01	VO.03.10.TI01
3	VO.03.10.SS02	VO.03.10.TI02
4	VO.03.10.SS03	VO.03.10.TI03
5	VO.03.10.SS04	VO.03.10.TI04
6	VO.03.10.SS05	VO.03.10.TI05
7	VO.03.10.SS06	VO.03.10.TI06
8	VO.03.10.SS07	VO.03.10.TI07
9	VO.03.10.SS08	VO.03.10.TI08
10	VO.03.10.SS09	VO.03.10.TI09
11	VO.03.10.SS10	VO.03.10.TI10
12	VO.03.10.SS11	VO.03.10.TI11

VO.03.9 3.5mm Small Dynamic Compression Plate

For Ulna and Radius, fixation with 3.5mm & 4.0mm Screws



Holes	Stainless Steel	Titanium
2	VO.03.9.SS01	VO.03.9.TI01
3	VO.03.9.SS02	VO.03.9.TI02
4	VO.03.9.SS03	VO.03.9.TI03
5	VO.03.9.SS04	VO.03.9.TI04
6	VO.03.9.SS05	VO.03.9.TI05
7	VO.03.9.SS06	VO.03.9.TI06
8	VO.03.9.SS07	VO.03.9.TI07
9	VO.03.9.SS08	VO.03.9.TI08
10	VO.03.9.SS09	VO.03.9.TI09
11	VO.03.9.SS10	VO.03.9.TI10
12	VO.03.9.SS11	VO.03.9.TI11
13	VO.03.9.SS12	VO.03.9.TI12
14	VO.03.9.SS13	VO.03.9.TI13

VO.03.11 3.5mm Small Limited Contact Dynamic Compression Plate

For Ulna and Radius, fixation with 3.5mm & 4.0mm Screws



Holes	Stainless Steel	Titanium
3	VO.03.11.SS01	VO.03.11.TI01
4	VO.03.11.SS02	VO.03.11.TI02
5	VO.03.11.SS03	VO.03.11.TI03
6	VO.03.11.SS04	VO.03.11.TI04
7	VO.03.11.SS05	VO.03.11.TI05
8	VO.03.11.SS06	VO.03.11.TI06
9	VO.03.11.SS07	VO.03.11.TI07
10	VO.03.11.SS08	VO.03.11.TI08
11	VO.03.11.SS09	VO.03.11.TI09
12	VO.03.11.SS10	VO.03.11.TI10
13	VO.03.11.SS11	VO.03.11.TI11
14	VO.03.11.SS12	VO.03.11.TI12

SECTION 3

SMALL FRAGMENT

Plates

VO.03.12 3.5mm Small T-Plate Right Angled with 3 Head Holes

For Distal Radius, Can be used for Left & Right Radius.
Fixation with 3.5mm & 4.0mm Screws



Holes	Stainless Steel	Titanium
3	VO.03.12.SS01	VO.03.12.TI01
4	VO.03.12.SS02	VO.03.12.TI02
5	VO.03.12.SS03	VO.03.12.TI03
6	VO.03.12.SS04	VO.03.12.TI04
7	VO.03.12.SS05	VO.03.12.TI05
8	VO.03.12.SS06	VO.03.12.TI06

VO.03.14, VO.03.15 3.5mm Small T-Plate, Oblique Angled

For Distal Radius, Can be used for Left & Right Radius.
Fixtation with 3.5mm & 4.0mm Screws

Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	VO.03.14.SS01L	VO.03.14.TI01L	VO.03.15.SS01R	VO.03.15.TI01R
4	VO.03.14.SS02L	VO.03.14.TI02L	VO.03.15.SS02R	VO.03.15.TI02R
5	VO.03.14.SS03L	VO.03.14.TI03L	VO.03.15.SS03R	VO.03.15.TI03R
6	VO.03.14.SS04L	VO.03.14.TI04L	VO.03.15.SS04R	VO.03.15.TI04R
7	VO.03.14.SS05L	VO.03.14.TI05L	VO.03.15.SS05R	VO.03.15.TI05R
8	VO.03.14.SS06L	VO.03.14.TI06L	VO.03.15.SS06R	VO.03.15.TI06R

VO.03.16 3.5mm Hook Plate



Holes	Stainless Steel	Titanium
4	VO.03.16.SS01	VO.03.16.TI01
5	VO.03.16.SS02	VO.03.16.TI02
6	VO.03.16.SS03	VO.03.16.TI03
7	VO.03.16.SS04	VO.03.16.TI04
8	VO.03.16.SS05	VO.03.16.TI05
10	VO.03.16.SS06	VO.03.16.TI06

VO.03.13 3.5mm Small T-Plate Right Angled with 4 Head Holes

For the Left & Right Radius.
Fixation with 3.5mm & 4.0mm Screws



Holes	Stainless Steel	Titanium
3	VO.03.13.SS01	VO.03.13.TI01
4	VO.03.13.SS02	VO.03.13.TI02
5	VO.03.13.SS03	VO.03.13.TI03
6	VO.03.13.SS04	VO.03.13.TI04
7	VO.03.13.SS05	VO.03.13.TI05
8	VO.03.13.SS06	VO.03.13.TI06



VO.03.17 3.5mm Pelvic Reconstruction Plate, Straight



Holes	Stainless Steel	Titanium
5	VO.03.17.SS01	VO.03.17.TI01
6	VO.03.17.SS02	VO.03.17.TI02
7	VO.03.17.SS03	VO.03.17.TI03
8	VO.03.17.SS04	VO.03.17.TI04
9	VO.03.17.SS05	VO.03.17.TI05
10	VO.03.17.SS06	VO.03.17.TI06
11	VO.03.17.SS07	VO.03.17.TI07
12	VO.03.17.SS08	VO.03.17.TI08
13	VO.03.17.SS09	VO.03.17.TI09
14	VO.03.17.SS10	VO.03.17.TI10
15	VO.03.17.SS11	VO.03.17.TI11
16	VO.03.17.SS12	VO.03.17.TI12
17	VO.03.17.SS13	VO.03.17.TI13
18	VO.03.17.SS14	VO.03.17.TI14
19	VO.03.17.SS15	VO.03.17.TI15
20	VO.03.17.SS16	VO.03.17.TI16
21	VO.03.17.SS17	VO.03.17.TI17
22	VO.03.17.SS18	VO.03.17.TI18

SECTION 3

SMALL FRAGMENT

Plates

VO.03.18


3.5mm Pelvic Reconstruction Plate, Curved



Holes	Stainless Steel	Titanium
5	VO.03.18.SS01	VO.03.18.TI01
6	VO.03.18.SS02	VO.03.18.TI02
7	VO.03.18.SS03	VO.03.18.TI03
8	VO.03.18.SS04	VO.03.18.TI04
9	VO.03.18.SS05	VO.03.18.TI05
10	VO.03.18.SS06	VO.03.18.TI06
11	VO.03.18.SS07	VO.03.18.TI07
12	VO.03.18.SS08	VO.03.18.TI08
13	VO.03.18.SS09	VO.03.18.TI09
14	VO.03.18.SS10	VO.03.18.TI10
16	VO.03.18.SS11	VO.03.18.TI11
18	VO.03.18.SS12	VO.03.18.TI12

VO.03.20

3.5mm Calcaneal Plate



Size (mm)	Stainless Steel	Titanium
60	VO.03.20.SS01	VO.03.20.TI01
70	VO.03.20.SS02	VO.03.20.TI02

VO.03.22

Cervical Spine H-Plate



Holes	Holes Spacing(mm)	Stainless Steel	Titanium
5	16	VO.03.22.SS01	VO.03.22.TI01
5	18	VO.03.22.SS02	VO.03.22.TI02
5	21	VO.03.22.SS03	VO.03.22.TI03
8	16+18	VO.03.22.SS04	VO.03.22.TI04
8	18+21	VO.03.22.SS05	VO.03.22.TI05
11	16+18+18	VO.03.22.SS06	VO.03.22.TI06

VO.03.19

3.5mm Semi Tubular Plate



Holes	Stainless Steel	Titanium
2	VO.03.19.SS01	VO.03.19.TI01
3	VO.03.19.SS02	VO.03.19.TI02
4	VO.03.19.SS03	VO.03.19.TI03
5	VO.03.19.SS04	VO.03.19.TI04
6	VO.03.19.SS05	VO.03.19.TI05
7	VO.03.19.SS06	VO.03.19.TI06
8	VO.03.19.SS07	VO.03.19.TI07
9	VO.03.19.SS08	VO.03.19.TI08
10	VO.03.19.SS09	VO.03.19.TI09
12	VO.03.19.SS10	VO.03.19.TI10

VO.03.21

3.5mm Cloverleaf Plate



Holes	Stainless Steel	Titanium
3	VO.03.21.SS01	VO.03.21.TI01
4	VO.03.21.SS02	VO.03.21.TI02
5	VO.03.21.SS03	VO.03.21.TI03
6	VO.03.21.SS04	VO.03.21.TI04
7	VO.03.21.SS05	VO.03.21.TI05
8	VO.03.21.SS06	VO.03.21.TI06
9	VO.03.21.SS07	VO.03.21.TI07
10	VO.03.21.SS08	VO.03.21.TI08
11	VO.03.21.SS09	VO.03.21.TI09
12	VO.03.21.SS10	VO.03.21.TI10
13	VO.03.21.SS11	VO.03.21.TI11
14	VO.03.21.SS12	VO.03.21.TI12

VO.03.23, VO.03.24
3.5mm Clavicle Hook Plate



Holes	Hook Depth (mm)	Left Direction		Right Direction	
		Stainless Steel	Titanium	Stainless Steel	Titanium
4	12 mm	VO.03.23.SS01L	VO.03.23.TI01L	VO.03.24.SS01R	VO.03.24.TI01R
4	15 mm	VO.03.23.SS02L	VO.03.23.TI02L	VO.03.24.SS02R	VO.03.24.TI02R
4	18 mm	VO.03.23.SS03L	VO.03.23.TI03L	VO.03.24.SS03R	VO.03.24.TI03R
5	12 mm	VO.03.23.SS04L	VO.03.23.TI04L	VO.03.24.SS04R	VO.03.24.TI04R
5	15 mm	VO.03.23.SS05L	VO.03.23.TI05L	VO.03.24.SS05R	VO.03.24.TI05R
5	18 mm	VO.03.23.SS06L	VO.03.23.TI06L	VO.03.24.SS06R	VO.03.24.TI06R
6	12 mm	VO.03.23.SS07L	VO.03.23.TI07L	VO.03.24.SS07R	VO.03.24.TI07R
6	15 mm	VO.03.23.SS08L	VO.03.23.TI08L	VO.03.24.SS08R	VO.03.24.TI08R
6	18 mm	VO.03.23.SS09L	VO.03.23.TI09L	VO.03.24.SS09R	VO.03.24.TI09R
7	12 mm	VO.03.23.SS10L	VO.03.23.TI10L	VO.03.24.SS10R	VO.03.24.TI10R
7	15 mm	VO.03.23.SS11L	VO.03.23.TI11L	VO.03.24.SS11R	VO.03.24.TI11R
7	18 mm	VO.03.23.SS12L	VO.03.23.TI12L	VO.03.24.SS12R	VO.03.24.TI12R

VO.03.25
3.5mm Y-Reconstruction Plate



Holes	Size	Stainless Steel	Titanium
6	Small	VO.03.25.SS01	VO.03.25.TI01
9	Medium	VO.03.25.SS02	VO.03.25.TI02
12	Large	VO.03.25.SS03	VO.03.25.TI03

Plates

VO.03.28, VO.03.29
Lecestre Plate

Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
5	VO.03.28.SS01L	VO.03.28.TI01L	VO.03.29.SS01R	VO.03.29.TI01R
7	VO.03.28.SS02L	VO.03.28.TI02L	VO.03.29.SS02R	VO.03.29.TI02R
9	VO.03.28.SS03L	VO.03.28.TI03L	VO.03.29.SS03R	VO.03.29.TI03R
10	VO.03.28.SS04L	VO.03.28.TI04L	VO.03.29.SS04R	VO.03.29.TI04R
11	VO.03.28.SS05L	VO.03.28.TI05L	VO.03.29.SS05R	VO.03.29.TI05R
13	VO.03.28.SS06L	VO.03.28.TI06L	VO.03.29.SS06R	VO.03.29.TI06R



VO.03.31
Proximal Humeral Plate Plate

Holes	Stainless Steel	Titanium
4	VO.03.31.SS01	VO.03.31.TI01
5	VO.03.31.SS02	VO.03.31.TI02
6	VO.03.31.SS03	VO.03.31.TI03
7	VO.03.31.SS04	VO.03.31.TI04
8	VO.03.31.SS05	VO.03.31.TI05
9	VO.03.31.SS06	VO.03.31.TI06
10	VO.03.31.SS07	VO.03.31.TI07



VO.03.33
Y-Fork Plate

Holes	Type	Stainless Steel	Titanium
4	Left	VO.03.33.SS01L	VO.03.33.TI01L
4	Right	VO.03.33.SS02R	VO.03.33.TI02R
5	Left	VO.03.33.SS03L	VO.03.33.TI03L
5	Right	VO.03.33.SS04R	VO.03.33.TI04R
6	Left	VO.03.33.SS05L	VO.03.33.TI05L
6	Right	VO.03.33.SS06R	VO.03.33.TI06R



SECTION 3

SMALL FRAGMENT

Instruments

VO.03.38 Drill Bit Quick Coupling End



Codes	Dia (mm)	Length (mm)
VO.03.38.01	2.5	110
VO.03.38.02	3.5	130

VO.03.39 Insert Drill Sleeve 3.5/2.5



VO.03.40 Tap



Codes	Dia (mm)
VO.03.40.01	3.5mm Cortical Screws
VO.03.40.02	4.0mm Cancellous Screws

VO.03.41 Double Drill Guide 3.5/2.5



VO.03.42 Universal Drill Guide 3.5mm



VO.03.43 Drill Guide 3.5, for Neutral and Loaded Position



VO.03.44 Counter Sink for 3.5/4.0mm Screws Quick Coupling



VO.03.45 Hexagonal Screw Driver Shaft 2.5mm Tip Quick Coupling



SECTION 3 SMALL FRAGMENT

Instruments

VO.03.46
Hexagonal Screw Driver with Holding
Sleeve 2.5mm Tip



VO.03.48
T-Handle Quick Coupling



VO.03.50
Sharp Hook



VO.03.52
Periosteal Elevator, Straight Edge 12mm



VO.03.47
Depth Gauge Measuring range up to 60mm



VO.03.49
Bending Iron (Pair)



VO.03.51
Periosteal Elevator, Round Edge - 6mm



VO.03.53
Hofmann Retractor, 8mm wide Short
Narrow Point, L-160mm



SECTION 3 SMALL FRAGMENT

Instruments


VO.03.56

Reduction Forceps, Serrated
Jaws Ratchet Lock, 140mm



VO.03.57

Reduction Forceps, Pointed,
Ratched Lock-140mm



VO.03.58

Small Fragment Instrument Set

VO.03.54

Hofmann Retractor, 15mm wide, L-



VO.03.55

Self Centering Bone Holding Forceps
Speed Lock, 190mm



Codes	Set Consisting of:	Units
VO.03.38.01	Drill Bit Quick Coupling End, Dia. 2.5mm x Length 110mm	2
VO.03.38.02	Drill Bit Quick Coupling End, Dia. 3.5mm x Length 130mm	2
VO.03.39	Insert Drill Sleeve 3.5/2.5	1
VO.03.40.01	Tap for Dia. 3.5mm Cortical Screws, Quick Coupling	2
VO.03.40.02	Tap for Dia. 4.0mm Cancellous Screws, Quick Coupling	2
VO.03.41	Double Drill Guide 3.5/2.5	1
VO.03.42	Universal Drill Guide 3.5	1
VO.03.43	Drill Guide 3.5, for Neutral and Loaded Position	1
VO.03.44	Counter Sink for 3.5/4.0mm Screws, Quick Coupling	1
VO.03.45	Hexagonal Screw Driver Shaft - 2.5mm Tip, Quick Coupling	1
VO.03.46	Hexagonal Screw Driver with Holding Sleeve-2.5mm Tip	1
VO.03.47	Depth Gauge, measuring range upto 60mm	1
VO.03.48	T-Handle Quick Coupling	1
VO.03.49	Bending Iron	1 Pair
VO.03.50	Sharp Hook	1
VO.03.51	Periosteal Elevator, Round Edge - 6mm	1
VO.03.52	Periosteal Elevator, Straight Edge - 12mm	1
VO.03.53	Hofmann Retractor, 8mm wide, Short Narrow Point, L-160mm	2
VO.03.54	Hofmann Retractor, 15mm wide, L-160mm	1
VO.03.55	Self Centering Bone Holding Forceps, Speed Lock, 190mm	2
VO.03.56	Reduction Forceps, Serrated Jaws, Speed lock, 140mm	1
VO.03.57	Reduction Forceps, Pointed, Ratchet Lock-140mm	1
VO.03.58.01	Aluminum Box Small Fragment Instrument Set	1

Instruments

VO.03.56

Reduction Forceps, Serrated
Jaws Ratchet Lock, 140mm

VO.03.54

Hofmann Retractor, 15mm wide, L-