

SECTION 5 (5.01-5.06)

Staples, Pins & Wires



**VO.05.5
Staples - Coventry**



Description	Codes
5mm Step	VO.05.1.SS01
10mm Step	VO.05.1.SS02
15mm Step	VO.05.1.SS03

**VO.05.6
Staples U - Type**



Dia (mm)	15 mm	20 mm	25 mm
1.5	VO.05.2.SS01	VO.05.2.SS05	VO.05.2.SS09
2.0	VO.05.2.SS02	VO.05.2.SS06	VO.05.2.SS10
2.5	VO.05.2.SS03	VO.05.2.SS07	VO.05.2.SS11
3.0	VO.05.2.SS04	VO.05.2.SS08	VO.05.2.SS12

Instruments for "U" Staples

**VO.05.3
U Staple - Driver**



**VO.05.5
U Staple - Inserter/Impactor**



**VO.05.7
U Staple - Extractor**



Instruments for Coventry Staples

**VO.05.4
Coventry Staple - Inserter**



**VO.05.6
Coventry Staple - Impactor**



**VO.05.8
Coventry Staple - Extractor**



VO.05.13
Suture Wire, 5 Meter - Stainless Steel



Codes	Size (mm)
VO.05.9.SS01	16 SWG (1.60mm)
VO.05.9.SS02	17 SWG (1.50mm)
VO.05.9.SS03	18 SWG (1.25mm)
VO.05.9.SS04	19 SWG (1.00mm)
VO.05.9.SS05	20 SWG (0.90mm)
VO.05.9.SS06	21 SWG (0.80mm)
VO.05.9.SS07	22 SWG (0.70mm)
VO.05.9.SS08	23 SWG (0.60mm)
VO.05.9.SS09	24 SWG (0.55mm)
VO.05.9.SS10	26 SWG (0.45mm)
VO.05.9.SS11	28 SWG (0.35mm)
VO.05.9.SS12	30 SWG (0.30mm)

VO.05.13
Suture Wire, 10 Meter - Stainless Steel



Codes	Size (mm)
VO.05.14.SS01	16 SWG (1.60mm)
VO.05.14.SS02	17 SWG (1.50mm)
VO.05.14.SS03	18 SWG (1.25mm)
VO.05.14.SS04	19 SWG (1.00mm)
VO.05.14.SS05	20 SWG (0.90mm)
VO.05.14.SS06	21 SWG (0.80mm)
VO.05.14.SS07	22 SWG (0.70mm)
VO.05.14.SS08	23 SWG (0.60mm)
VO.05.14.SS09	24 SWG (0.55mm)
VO.05.14.SS10	26 SWG (0.45mm)
VO.05.14.SS11	28 SWG (0.35mm)
VO.05.14.SS12	30 SWG (0.30mm)

VO.05.10
Kirschner Wire with Trocar Tip Both End



Diameter (mm)	Length				
	70 mm	100 mm	150 mm	225 mm	300 mm
0.8	VO.05.10.SS01	VO.05.10.SS15	VO.05.10.SS29	VO.05.10.SS43	VO.05.10.SS57
0.9	VO.05.10.SS02	VO.05.10.SS16	VO.05.10.SS30	VO.05.10.SS44	VO.05.10.SS58
1.0	VO.05.10.SS03	VO.05.10.SS17	VO.05.10.SS31	VO.05.10.SS45	VO.05.10.SS59
1.1	VO.05.10.SS04	VO.05.10.SS18	VO.05.10.SS32	VO.05.10.SS46	VO.05.10.SS60
1.2	VO.05.10.SS05	VO.05.10.SS19	VO.05.10.SS33	VO.05.10.SS47	VO.05.10.SS61
1.4	VO.05.10.SS06	VO.05.10.SS20	VO.05.10.SS34	VO.05.10.SS48	VO.05.10.SS62
1.5	VO.05.10.SS07	VO.05.10.SS21	VO.05.10.SS35	VO.05.10.SS49	VO.05.10.SS63
1.6	VO.05.10.SS08	VO.05.10.SS22	VO.05.10.SS36	VO.05.10.SS50	VO.05.10.SS64
1.8	VO.05.10.SS09	VO.05.10.SS23	VO.05.10.SS37	VO.05.10.SS51	VO.05.10.SS65
2.0	VO.05.10.SS10	VO.05.10.SS24	VO.05.10.SS38	VO.05.10.SS52	VO.05.10.SS66
2.2	VO.05.10.SS11	VO.05.10.SS25	VO.05.10.SS39	VO.05.10.SS53	VO.05.10.SS67
2.3	VO.05.10.SS12	VO.05.10.SS26	VO.05.10.SS40	VO.05.10.SS54	VO.05.10.SS68
2.5	VO.05.10.SS13	VO.05.10.SS27	VO.05.10.SS41	VO.05.10.SS55	VO.05.10.SS69
3.0	VO.05.10.SS14	VO.05.10.SS28	VO.05.10.SS42	VO.05.10.SS56	VO.05.10.SS70

VO.05.11 Kirschner Wire Fully Threaded



Codes	Dia (mm)	Length (mm)
VO.05.11.SS01	1.5	225
VO.05.11.SS02	2.0	225
VO.05.11.SS03	2.5	225
VO.05.11.SS04	1.5	150
VO.05.11.SS05	2.0	150
VO.05.11.SS06	2.5	150

VO.05.12 Steinman Pins



Length (mm)	Diameters				
	3 mm	3.5 mm	4 mm	4.5 mm	5 mm
150	VO.05.12.SS01	VO.05.12.SS05	VO.05.12.SS09	VO.05.12.SS13	VO.05.12.SS17
175	VO.05.12.SS02	VO.05.12.SS06	VO.05.12.SS10	VO.05.12.SS14	VO.05.12.SS18
200	VO.05.12.SS03	VO.05.12.SS07	VO.05.12.SS11	VO.05.12.SS15	VO.05.12.SS19
225	VO.05.12.SS04	VO.05.12.SS08	VO.05.12.SS12	VO.05.12.SS16	VO.05.12.SS20

VO.05.13 Centrally Threaded Steinman Pins



Codes	Dia (mm)	Length (mm)
VO.05.13.SS01	4.0	225
VO.05.13.SS02	4.5	225
VO.05.13.SS03	5.0	225
VO.05.13.SS04	6.0	225

VO.05. Circlage Wire With Loop



Size (mm)	Length 280mm	Length 600mm
17SWG(1.50mm)	VO.05.14.SS01	VO.05.14.SS06
18SWG(1.25mm)	VO.05.14.SS02	VO.05.14.SS07
19SWG(1.00mm)	VO.05.14.SS03	VO.05.14.SS08
20SWG(0.90mm)	VO.05.14.SS04	VO.05.14.SS09
21SWG(0.80mm)	VO.05.14.SS05	VO.05.14.SS10

VO.05.14 Austin Moore Pins, 3.0mm Dia with 2 Nuts



Codes	Length (mm)
VO.05.15.SS01	50
VO.05.15.SS02	55
VO.05.15.SS03	60
VO.05.15.SS04	65
VO.05.15.SS05	70
VO.05.15.SS06	75
VO.05.15.SS07	80
VO.05.15.SS08	85
VO.05.15.SS09	90
VO.05.15.SS10	95
VO.05.15.SS11	100
VO.05.15.SS12	105
VO.05.15.SS13	110
VO.05.15.SS14	115
VO.05.15.SS15	120
VO.05.15.SS16	125
VO.05.15.SS17	130
VO.05.15.SS18	135
VO.05.15.SS19	140
VO.05.15.SS20	145
VO.05.15.SS21	150

VO.05.15 Knowles Pin



Codes	Dia (mm)	Length (mm)
VO.05.16.SS01	3	65
VO.05.16.SS02	3	70
VO.05.16.SS03	3	75
VO.05.16.SS04	3	80
VO.05.16.SS05	3	85
VO.05.16.SS06	3	90
VO.05.16.SS07	3	95
VO.05.16.SS08	3	100
VO.05.16.SS09	4	65
VO.05.16.SS10	4	70
VO.05.16.SS11	4	75
VO.05.16.SS12	4	80
VO.05.16.SS13	4	85
VO.05.16.SS14	4	90
VO.05.16.SS15	4	95
VO.05.16.SS16	4	100

VO.05.16 Drummer Wire



Codes	Size
VO.05.17.SS01	18 SWG (1.25mm)
VO.05.17.SS02	19 SWG (1.00mm)
VO.05.17.SS03	20 SWG (0.90mm)

VO.05.18 Wire Implant Set

Code	Set Consisting of:	Units
VO.05.10.SS46	K.Wire Dia 1.0mm, Length 225mm	15
VO.05.10.SS49	K.Wire Dia 1.5mm, Length 225mm	15
VO.05.10.SS52	K.Wire Dia 2.0mm, Length 225mm	15
VO.05.10.SS55	K.Wire Dia 2.5mm, Length 225mm	15
VO.05.10.SS56	K.Wire Dia 3.0mm, Length 225mm	15
VO.05.9.SS03	S.S. Wire : 18 SWG (2 each)	2
VO.05.9.SS05	S.S. Wire : 20 SWG (2 each)	2
VO.05.9.SS07	S.S. Wire : 22 SWG (2 each)	2
VO.05.12.SS08	Steinmann Pin Dia 3.5, Length 225mm (5 each)	5
VO.05.12.SS12	Steinmann Pin Dia 4.0, Length 225mm (5 each)	5
VO.05.12.SS16	Steinmann Pin Dia 4.5, Length 225mm (5 each)	5
VO.05.18.01	Aluminum Box for Wire Implant Set	1

Instruments

VO.05.19
Wire Tightener & Twister



VO.05.20
Wire Tightner cum Twister & Cutter



VO.05.21
K. Wire Introducer



VO.05.22
K. Wire Bender



VO.05.23
Triple Drill guide 2.0mm



VO.05.24
Bone Awl



VO.05.25
K. Wire Tractor



VO.05.26
Wire Passer



Codes	Size
VO.05.25.01	Small
VO.05.25.02	Medium
VO.05.25.03	Large

Codes	Size
VO.05.26.01	Small
VO.05.26.02	Medium
VO.05.26.03	Large

VO.05.27
Plier cum Wire Bender & Wire Cutter

VO.05.28
Circlage Wire Tightner



Codes	Size (inch)
VO.05.27.01	9
VO.05.27.02	11

VO.05.29
T- Handle with Chuck & Key



VO.05.31
Steinmann Pin Holder



VO.05.30
Steinmann Pin Cutter / Plate Cutter
(Capacity 5mm)



Codes	Size
VO.05.31.01	Small
VO.05.31.02	Medium
VO.05.31.03	Large

VO.05.32
Wire Instrument Set

Codes	Set Consisting of:	Units
VO.05.19	Wire Tightener Cum Twister	1
VO.05.27.01	Plier cum Wire Bender cum Wire Cutter 9"	1



Innovating Orthopedic Care

SECTION 5 STAPLES, PINS & WIRES

VO.05.20	K.Wire Bender	1
VO.05.26.02	Wire Passer - Medium	2
VO.05.21	K.Wire Introducer	1
VO.05.29	T- Handle with Chuck (Steinmann Pin Introducer)	1
VO.05.24	Bone Awl with Eye	1
VO.05.32.01	Aluminum Box for Wire Instrument Set	1